

BCL2 / Bcl-2 Rabbit anti-Human Polyclonal (pSer70) (FITC) Antibody - LS-C252480 - LSBio	
CatalogID:	LS-C252480
Target:	B-cell CLL/lymphoma 2 (BCL2)
Synonyms:	BCL2 Antibody, Apoptosis regulator Bcl-2 Antibody, B-cell CLL/lymphoma 2 Antibody, B-cell lymphoma protein 2 Antibody, Bcl-2 Antibody, PPP1R50 Antibody
Family / Subfamily:	Apoptosis
Host	BCL2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	BCL2 / Bcl-2 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	BCL2 / Bcl-2 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser70 of human bcl.
Specificity:	Human BCL2 / Bcl-2
Epitope:	pSer70
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014 © 2014 LifeSpan BioSciences	