

BMP7 / BMP-7 Rabbit anti-Human Polyclonal (aa279-309) (PE) Antibody - LS-C252617 - LSBio	
CatalogID:	LS-C252617
Target:	bone morphogenetic protein 7 (BMP7)
Synonyms:	BMP7 Antibody, Bone morphogenetic protein 7 Antibody, Eptotermin alfa Antibody, OP-1 Antibody, OP1 Antibody, Osteogenic protein 1 Antibody, BMP-7 Antibody
Family / Subfamily:	TGF beta / not assigned-TGF beta
Host	BMP7 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	BMP7 / BMP-7 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	BMP7 / BMP-7 antibody was raised against kLH-conjugated synthetic peptide (aa279-309) from the C-terminus of human BMP7.
Specificity:	Human BMP7 / BMP-7
Epitope:	aa279-309
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	