

NUAK1 / ARK5 Rabbit anti-Human Polyclonal (aa4-33) (PE) Antibody - LS-C252883 - LSBio	
CatalogID:	LS-C252883
Target:	NUAK family, SNF1-like kinase, 1 (NUAK1)
Synonyms:	NUAK1 Antibody, KIAA0537 Antibody, NUAKE family SNF1-like kinase 1 Antibody, OMPHK1 Antibody, AMPK-related protein kinase 5 Antibody, ARK5 Antibody, Omphalocele kinase 1 Antibody
Family / Subfamily:	Protein Kinase / CAMKL
Host	NUAK1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	NUAK1 / ARK5 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NUAK1 / ARK5 antibody was raised against KLH-conjugated synthetic peptide (aa4-33) from the N-terminal region of human ARK5.
Specificity:	Human NUAKE1 / ARK5.
Epitope:	aa4-33
Reactivity:	Human, Mouse, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	