

BPGM Rabbit anti-Human Polyclonal (aa229-259) (FITC) Antibody - LS-C252974 - LSBio	
CatalogID:	LS-C252974
Target:	2,3-bisphosphoglycerate mutase (BPGM)
Synonyms:	BPGM Antibody, 2,3-diphosphoglycerate mutase Antibody, 2,3-bisphosphoglycerate mutase Antibody, Bisphosphoglycerate mutase Antibody, BPG-dependent PGAM Antibody, DPGM Antibody
Host	BPGM antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	BPGM antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	BPGM antibody was raised against kLH-conjugated synthetic peptide (aa229-259) from the C-terminus of human BPGM.
Specificity:	Human BPGM
Epitope:	aa229-259
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	