

TSPEAR Rabbit anti-Human Polyclonal (aa299-329) (PE) Antibody - LS-C253156 - LSBio	
CatalogID:	LS-C253156
Target:	Thrombospondin-type laminin G domain and EAR repeats (TSPEAR)
Synonyms:	TSPEAR Antibody
Host	TSPEAR antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	TSPEAR antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TSPEAR antibody was raised against kLH-conjugated synthetic peptide (aa299-329) from internal region of human C21orf29.
Specificity:	Human TSPEAR
Epitope:	aa299-329
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	