

CGRP / Calcitonin Rabbit anti-Human Polyclonal (aa63-91) (FITC) Antibody - LS-C253298 - LSBio	
CatalogID:	LS-C253298
Target:	calcitonin-related polypeptide alpha (CALCA)
Synonyms:	CALCA Antibody, Alpha-type CGRP Antibody, CGRP Antibody, CGRP-I Antibody, CGRP1 Antibody, CT Antibody, Calcitonin Antibody, KC Antibody, CALC1 Antibody, Calcitonin 1 Antibody, Katalcalcin Antibody
Family / Subfamily:	Hormone / not assigned-Hormone
Host	CALCA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CGRP / Calcitonin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CGRP / Calcitonin antibody was raised against kLH-conjugated synthetic peptide (aa63-91) from internal region of human CALCA.
Specificity:	Human CGRP / Calcitonin
Epitope:	aa63-91
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	