

CASP3 / Caspase 3 Rabbit anti-Human Polyclonal (Asp175) (FITC) Antibody - LS-C253407 - LSBio	
CatalogID:	LS-C253407
Target:	caspase 3, apoptosis-related cysteine peptidase (CASP3)
Synonyms:	CASP3 Antibody, Apopain Antibody, CASP-3 Antibody, Caspase-3 Antibody, CPP32B Antibody, CPP32 Antibody, Cysteine protease CPP32 Antibody, SCA-1 Antibody, SREBP cleavage activity 1 Antibody, Yama Antibody, Procaspase3 Antibody, Caspase 3 Antibody, CPP-32 Antibody, PARP cleavage protease Antibody, Protein Yama Antibody
Family / Subfamily:	Protease / Cysteine C14
Host	CASP3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CASP3 / Caspase 3 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CASP3 / Caspase 3 antibody was raised against kLH-conjugated synthetic phosphopeptide (aa155-184) surrounding Asp175 of human CASP3.
Specificity:	Human CASP3 / Caspase 3
Epitope:	Asp175
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	