

CASP3 / Caspase 3 Rabbit anti-Mouse Polyclonal (pSer26) (FITC) Antibody - LS-C253418 - LSBio	
CatalogID:	LS-C253418
Target:	caspase 3, apoptosis-related cysteine peptidase (CASP3)
Synonyms:	CASP3 Antibody, Apopain Antibody, CASP-3 Antibody, Caspase-3 Antibody, CPP32B Antibody, CPP32 Antibody, Cysteine protease CPP32 Antibody, SCA-1 Antibody, SREBP cleavage activity 1 Antibody, Yama Antibody, Procaspase3 Antibody, Caspase 3 Antibody, CPP-32 Antibody, PARP cleavage protease Antibody, Protein Yama Antibody
Family / Subfamily:	Protease / Cysteine C14
Host	CASP3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CASP3 / Caspase 3 antibody was raised against Mouse
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CASP3 / Caspase 3 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser26 of mouse CASP3.
Specificity:	Mouse CASP3 / Caspase 3
Epitope:	pSer26
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014 © 2014 LifeSpan BioSciences	