

**CASP3 / Caspase 3 Rabbit anti-Human Polyclonal (aa155-185) (FITC) Antibody - LS-C253425 - LSBio**

<b>CatalogID:</b>	LS-C253425
<b>Target:</b>	caspase 3, apoptosis-related cysteine peptidase (CASP3)
<b>Synonyms:</b>	CASP3 Antibody, Apopain Antibody, CASP-3 Antibody, Caspase-3 Antibody, CPP32B Antibody, CPP32 Antibody, Cysteine protease CPP32 Antibody, SCA-1 Antibody, SREBP cleavage activity 1 Antibody, Yama Antibody, Procaspase3 Antibody, Caspase 3 Antibody, CPP-32 Antibody, PARP cleavage protease Antibody, Protein Yama Antibody
<b>Family / Subfamily:</b>	Protease / Cysteine C14
<b>Host</b>	CASP3 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CASP3 / Caspase 3 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CASP3 / Caspase 3 antibody was raised against kLH-conjugated synthetic peptide (aa155-185) from the Asp175 region of human CASP3.
<b>Specificity:</b>	Human CASP3 / Caspase 3
<b>Epitope:</b>	aa155-185
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences