

CASP6 / Caspase 6 Rabbit anti-Human Polyclonal (aa16-45) (PE) Antibody - LS-C253441 - LSBio	
CatalogID:	LS-C253441
Target:	caspase 6, apoptosis-related cysteine peptidase (CASP6)
Synonyms:	CASP6 Antibody, Apoptotic protease Mch-2 Antibody, Caspase 6 Antibody, CASP-6 Antibody, Caspase-6 Antibody, Mch-2 Antibody, MCH2 Antibody
Family / Subfamily:	Protease / Cysteine C14
Host	CASP6 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CASP6 / Caspase 6 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CASP6 / Caspase 6 antibody was raised against kLH-conjugated synthetic peptide (aa16-45) from the N-terminal region of human CASP6.
Specificity:	Human CASP6 / Caspase 6
Epitope:	aa16-45
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	