

CCNI2 Rabbit anti-Human Polyclonal (aa99-129) (PE) Antibody - LS-C253768 - LSBio

CatalogID:	LS-C253768
Target:	cyclin I family, member 2 (CCNI2)
Synonyms:	CCNI2 Antibody, Cyclin I family, member 2 Antibody, Cyclin-I2 Antibody
Host	CCNI2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CCNI2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CCNI2 antibody was raised against kLH-conjugated synthetic peptide (aa99-129) from internal region of human CCNI2.
Specificity:	Human CCNI2
Epitope:	aa99-129
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences