

CCR5 Rabbit anti-Human Polyclonal (aa151-178) (FITC) Antibody - LS-C253814 - LSBio

CatalogID:	LS-C253814
Target:	chemokine (C-C motif) receptor 5 (gene/pseudogene) (CCR5)
Synonyms:	CCR5 Antibody, C-C chemokine receptor 5 Antibody, CD195 antigen Antibody, CD195 Antibody, C-C CKR-5 Antibody, Chemokine c-c motif receptor 5 Antibody, Chemr13 Antibody, CMKBR5 Antibody, CKR5 Antibody, Cc chemokine receptor 5 Antibody, CCCR5 Antibody, CCR-5 Antibody, Ccr5 receptor Antibody, Chemokine receptor CCR5 Antibody, HIV-1 fusion coreceptor Antibody, IDDM22 Antibody, C-C chemokine receptor type 5 Antibody, CC-CKR-5 Antibody, CKR-5 Antibody
Family / Subfamily:	GPCR / Chemokine
Host	CCR5 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CCR5 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CCR5 antibody was raised against kLH-conjugated synthetic peptide (aa151-178) from internal region of human CCR5.
Specificity:	Human CCR5.
Epitope:	aa151-178
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences