

CD19 Rabbit anti-Human Polyclonal (aa400-429) (Azide-free) (HRP) Antibody - LS-C253895 - LSBio

CatalogID:	LS-C253895
Target:	CD19 molecule
Synonyms:	CD19 Antibody, B-lymphocyte antigen CD19 Antibody, CD19 molecule Antibody, Differentiation antigen CD19 Antibody, T-cell surface antigen Leu-12 Antibody, B4 Antibody, CD19 antigen Antibody, CVID3 Antibody
Host	CD19 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	CD19 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD19 antibody was raised against kLH-conjugated synthetic peptide (aa400-429) from the C-terminus of human CD19.
Specificity:	Human CD19.
Epitope:	aa400-429
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, azide-free.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences