

CD1A Rabbit anti-Human Polyclonal (aa74-104) (FITC) Antibody - LS-C253899 - LSBio

CatalogID:	LS-C253899
Target:	CD1a molecule
Synonyms:	CD1A Antibody, CD1 Antibody, CD1A antigen, a polypeptide Antibody, CD1a molecule Antibody, FCB6 Antibody, HTA1 Antibody, HTa1 thymocyte antigen Antibody, R4 Antibody, T6 Antibody, CD1a antigen Antibody, Cluster of differentiation 1 A Antibody
Host	CD1A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD1A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD1A antibody was raised against kLH-conjugated synthetic peptide (aa74-104) from internal region of human CD1A.
Specificity:	Human CD1A.
Epitope:	aa74-104
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences