

CD5 Rabbit anti-Human Polyclonal (aa202-232) (FITC) Antibody - LS-C254078 - LSBio

CatalogID:	LS-C254078
Target:	CD5 molecule
Synonyms:	CD5 Antibody, CD5 antigen Antibody, CD5 antigen (p56-62) Antibody, CD5 molecule Antibody, Lymphocyte antigen T1/Leu-1 Antibody, LEU1 Antibody, T1 Antibody
Host	CD5 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD5 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD5 antibody was raised against kLH-conjugated synthetic peptide (aa202-232) from internal region of human CD5.
Specificity:	Human CD5.
Epitope:	aa202-232
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences