

CD7 Rabbit anti-Human Polyclonal (aa62-92) (FITC) Antibody - LS-C254090 - LSBio

CatalogID:	LS-C254090
Target:	CD7 molecule
Synonyms:	CD7 Antibody, CD7 antigen Antibody, CD7 antigen (p41) Antibody, CD7 molecule Antibody, LEU-9 Antibody, T-cell leukemia antigen Antibody, T-cell antigen CD7 Antibody, T-cell surface antigen Leu-9 Antibody, TP41 Antibody, gp40 Antibody, p41 protein Antibody, Tp40 Antibody
Host	CD7 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD7 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD7 antibody was raised against kLH-conjugated synthetic peptide (aa62-92) from internal region of human CD7.
Specificity:	Human CD7.
Epitope:	aa62-92
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences