

**CD71 / Transferrin Receptor Rabbit anti-Human Polyclonal (aa656-685) (PE) Antibody - LS-C254112 - LSBio**

<b>CatalogID:</b>	LS-C254112
<b>Target:</b>	transferrin receptor (TFRC)
<b>Synonyms:</b>	TFRC Antibody, CD71 Antibody, TFR1 Antibody, TRFR Antibody, TFR Antibody, TR Antibody, CD71 antigen Antibody, p90 Antibody, T9 Antibody, Transferrin receptor Antibody, Transferrin receptor protein 1 Antibody
<b>Family / Subfamily:</b>	Protease / Metallopeptidase M28B
<b>Host</b>	TFRC antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD71 / Transferrin Receptor antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CD71 / Transferrin Receptor antibody was raised against kLH-conjugated synthetic peptide (aa656-685) from the C-terminus of human CD71.
<b>Specificity:</b>	Human CD71 / Transferrin Receptor
<b>Epitope:</b>	aa656-685
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences