

CD8A / CD8 Alpha Rabbit anti-Human Polyclonal (aa158-188) (PE) Antibody - LS-C254145 - LSBio	
CatalogID:	LS-C254145
Target:	CD8a molecule
Synonyms:	CD8A Antibody, CD8a molecule Antibody, CD8a antigen Antibody, OKT8 T-cell antigen Antibody, Leu2 T-lymphocyte antigen Antibody, T8 T-cell antigen Antibody, T cell co-receptor Antibody, p32 Antibody, CD8 Antibody, Leu2 Antibody, T-cell antigen Leu2 Antibody
Host	CD8A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD8A / CD8 Alpha antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD8A / CD8 Alpha antibody was raised against kLH-conjugated synthetic peptide (aa158-188) from the C-terminus of human CD8A.
Specificity:	Human CD8A / CD8 Alpha
Epitope:	aa158-188
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014 © 2014 LifeSpan BioSciences	