

CD8A / CD8 Alpha Rabbit anti-Human Polyclonal (aa66-96) (FITC) Antibody - LS-C254151 - LSBio	
<b>CatalogID:</b>	LS-C254151
<b>Target:</b>	CD8a molecule
<b>Synonyms:</b>	CD8A Antibody, CD8a molecule Antibody, CD8a antigen Antibody, OKT8 T-cell antigen Antibody, Leu2 T-lymphocyte antigen Antibody, T8 T-cell antigen Antibody, T cell co-receptor Antibody, p32 Antibody, CD8 Antibody, Leu2 Antibody, T-cell antigen Leu2 Antibody
<b>Host</b>	CD8A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD8A / CD8 Alpha antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CD8A / CD8 Alpha antibody was raised against kLH-conjugated synthetic peptide (aa66-96) from the N-terminal region of human CD8A.
<b>Specificity:</b>	Human CD8A / CD8 Alpha
<b>Epitope:</b>	aa66-96
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b> Created on 9/29/2014 © 2014 LifeSpan BioSciences	