

CD9 Rabbit anti-Human Polyclonal (aa29-58) (FITC) Antibody - LS-C254167 - LSBio

CatalogID:	LS-C254167
Target:	CD9 molecule
Synonyms:	CD9 Antibody, 5H9 antigen Antibody, BA2 Antibody, BA-2/p24 antigen Antibody, CD9 antigen Antibody, CD9 antigen (p24) Antibody, DRAP-27 Antibody, MIC3 Antibody, MRP-1 Antibody, Leukocyte antigen MIC3 Antibody, Tetraspanin-29 Antibody, TSPAN29 Antibody, BTCC-1 Antibody, CD9 molecule Antibody, Motility related protein-1 Antibody, Motility-related protein Antibody, TSPAN-29 Antibody
Family / Subfamily:	Tetraspan / not assigned-Tetraspan
Host	CD9 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD9 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD9 antibody was raised against kLH-conjugated synthetic peptide (aa29-58) from the N-terminal region of human CD9.
Specificity:	Human CD9
Epitope:	aa29-58
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences