

CDC25A Rabbit anti-Human Polyclonal (Ser78) (FITC) Antibody - LS-C254227 - LSBio	
<b>CatalogID:</b>	LS-C254227
<b>Target:</b>	cell division cycle 25A (CDC25A)
<b>Synonyms:</b>	CDC25A Antibody, CDC25A2-CAG isoform Antibody, Cell division cycle 25A Antibody, CDC25A phosphatase Antibody, CDC25Hu1 Antibody, M-phase inducer phosphatase 1 Antibody, CDC25A2 Antibody
<b>Family / Subfamily:</b>	Protein Phosphatase / CDC25
<b>Host</b>	CDC25A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CDC25A antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CDC25A antibody was raised against kLH-conjugated synthetic peptide around Ser78 of human CDC25A.
<b>Specificity:</b>	Human CDC25A
<b>Epitope:</b>	Ser78
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	