

CDC25B Rabbit anti-Human Polyclonal (Ser353) (FITC) Antibody - LS-C254245 - LSBio

CatalogID:	LS-C254245
Target:	cell division cycle 25B (CDC25B)
Synonyms:	CDC25B Antibody, CDC25B phosphatase Antibody, Cell division cycle 25B Antibody, M-phase inducer phosphatase 2 Antibody, CDC25HU2 Antibody
Family / Subfamily:	Protein Phosphatase / CDC25
Host	CDC25B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CDC25B antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CDC25B antibody was raised against kLH-conjugated synthetic peptide around Ser353 of human CDC25B.
Specificity:	Human CDC25B.
Epitope:	Ser353
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences