

CDC2 / CDK2 Rabbit anti-Mouse Polyclonal (aa226-254) (FITC) Antibody - LS-C254423 - LSBio	
<b>CatalogID:</b>	LS-C254423
<b>Target:</b>	cyclin-dependent kinase 2 (CDK2)
<b>Synonyms:</b>	CDK2 Antibody, CDC2 Antibody, Cell division kinase 2 Antibody, Cell division protein kinase 2 Antibody, CDC2-related protein kinase Antibody, Cell division kinase 2 Antibody, p33 protein kinase Antibody, p33(CDK2) Antibody, CDKN2 Antibody, Cyclin-dependent kinase 2 Antibody
<b>Family / Subfamily:</b>	Protein Kinase / CDC2/CDK
<b>Host</b>	CDK2 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CDC2 / CDK2 antibody was raised against Mouse
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CDC2 / CDK2 antibody was raised against KLH-conjugated synthetic peptide (aa226-254) from the C-terminus of Mouse Cdk2.
<b>Specificity:</b>	Mouse CDC2 / CDK2
<b>Epitope:</b>	aa226-254
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	