

CDX2 Rabbit anti-Human Polyclonal (Ser283) (PE) Antibody - LS-C254544 - LSBio

CatalogID:	LS-C254544
Target:	caudal type homeobox 2 (CDX2)
Synonyms:	CDX2 Antibody, Caudal-type homeobox protein 2 Antibody, CDX3 Antibody, Homeobox protein CDX-2 Antibody, Caudal type homeobox 2 Antibody, CDX-3 Antibody
Host	CDX2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CDX2 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CDX2 antibody was raised against kLH-conjugated synthetic phosphopeptide (aa268-297) around Ser283 of human CDX2.
Specificity:	Human CDX2.
Epitope:	Ser283
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences