

CEP70 Rabbit anti-Human Polyclonal (aa316-346) (PE) Antibody - LS-C254643 - LSBio

CatalogID:	LS-C254643
Target:	centrosomal protein 70kDa (CEP70)
Synonyms:	CEP70 Antibody, Centrosomal protein 70kDa Antibody, Centrosomal protein of 70 kDa Antibody, p10-binding protein Antibody, BITE Antibody
Host	CEP70 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CEP70 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CEP70 antibody was raised against KLH-conjugated synthetic peptide (aa316-346) from internal region of human CEP70.
Specificity:	Human CEP70
Epitope:	aa316-346
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences