

CLTC / Clathrin Heavy Chain Rabbit anti-Human Polyclonal (aa1024-1054) (PE) Antibody - LS-C255096 - LSBio

CatalogID:	LS-C255096
Target:	clathrin, heavy chain (Hc) (CLTC)
Synonyms:	CLTC Antibody, Clathrin, heavy chain Antibody, CLH17 Antibody, CLTCL2 Antibody, CHC Antibody, CHC17 Antibody, Clathrin heavy chain 1 Antibody, CLH-17 Antibody, KIAA0034 Antibody, Hc Antibody, Clathrin Heavy Chain Antibody, Clathrin, heavy chain (Hc) Antibody
Host	CLTC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CLTC / Clathrin Heavy Chain antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CLTC / Clathrin Heavy Chain antibody was raised against kLH-conjugated synthetic peptide (aa1024-1054) from internal region of human CLTC.
Specificity:	Human CLTC / Clathrin Heavy Chain
Epitope:	aa1024-1054
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences