

COBL Rabbit anti-Human Polyclonal (aa191-219) (FITC) Antibody - LS-C255176 - LSBio

CatalogID:	LS-C255176
Target:	cordon-bleu WH2 repeat protein (COBL)
Synonyms:	COBL Antibody, Cordon-bleu homolog (mouse) Antibody, Cordon-bleu WH2 repeat protein Antibody, KIAA0633 Antibody, Protein cordon-bleu Antibody
Host	COBL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	COBL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	COBL antibody was raised against KLH-conjugated synthetic peptide (aa191-219) from the N-terminal region of human COBL.
Specificity:	Human COBL
Epitope:	aa191-219
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences