

C16orf13 Rabbit anti-Human Polyclonal (aa1-30) (FITC) Antibody - LS-C255299 - LSBio

CatalogID:	LS-C255299
Target:	chromosome 16 open reading frame 13 (C16orf13)
Synonyms:	C16orf13 Antibody, JFP2 Antibody, MGC13114 Antibody, UPF0585 protein C16orf13 Antibody
Host	C16orf13 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	C16orf13 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	C16orf13 antibody was raised against kLH-conjugated synthetic peptide (aa1-30) from the N-terminal region of human CP013.
Specificity:	Human C16orf13.
Epitope:	aa1-30
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences