

CRYM Rabbit anti-Mouse Polyclonal (aa154-183) (FITC) Antibody - LS-C255493 - LSBio

CatalogID:	LS-C255493
Target:	crystallin, mu (CRYM)
Synonyms:	CRYM Antibody, Crystallin, mu Antibody, Ketimine reductase Antibody, Mu-crystallin homolog Antibody, THBP Antibody, DFNA40 Antibody
Host	CRYM antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CRYM antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	CRYM antibody was raised against kLH-conjugated synthetic peptide (aa154-183) from internal region of human CRYM.
Specificity:	Mouse CRYM.
Epitope:	aa154-183
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences