

CTNNB1 / Beta Catenin Rabbit anti-Human Polyclonal (aa653-681) (PE) Antibody - LS-C255650 - LSBio	
CatalogID:	LS-C255650
Target:	catenin (cadherin-associated protein), beta 1, 88kDa (CTNNB1)
Synonyms:	CTNNB1 Antibody, Beta-catenin Antibody, Catenin beta-1 Antibody, CTNNB Antibody, Beta Catenin Antibody
Family / Subfamily:	Catenin / not assigned-Catenin
Host	CTNNB1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CTNNB1 / Beta Catenin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CTNNB1 / Beta Catenin antibody was raised against kLH-conjugated synthetic peptide (aa653-681) from the C-terminus of human CTNNB1.
Specificity:	Human CTNNB1 / Beta Catenin
Epitope:	aa653-681
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	