

| DMGDH Rabbit anti-Human Polyclonal (aa835-864) (FITC) Antibody - LS-C256148 - LSBio | |
|---|---|
| CatalogID: | LS-C256148 |
| Target: | dimethylglycine dehydrogenase (DMGDH) |
| Synonyms: | DMGDH Antibody, Dimethylglycine dehydrogenase Antibody, DMGDHD Antibody, ME2GLYDH Antibody |
| Host | DMGDH antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | DMGDH antibody was raised against Human |
| Antigen Type: | Synthetic peptide |
| Immunogen: | DMGDH antibody was raised against kLH-conjugated synthetic peptide (aa835-864) from the C-terminus of human DMGDH. |
| Specificity: | Human DMGDH |
| Epitope: | aa835-864 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/29/2014 © 2014 LifeSpan BioSciences | |