

E2F1 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C256426 - LSBio

CatalogID:	LS-C256426
Target:	E2F transcription factor 1 (E2F1)
Synonyms:	E2F1 Antibody, E2F transcription factor 1 Antibody, PRB-binding protein E2F-1 Antibody, RBAP-1 Antibody, RBBP3 Antibody, RBBP-3 Antibody, Transcription factor E2F1 Antibody, RBAP1 Antibody, E2F-1 Antibody, PBR3 Antibody
Host	E2F1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	E2F1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	E2F1 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around H357 of human E2F1.
Specificity:	Human E2F1.
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences