

EBPL Rabbit anti-Human Polyclonal (aa184-214) (PE) Antibody - LS-C256462 - LSBio

CatalogID:	LS-C256462
Target:	Emopamil Binding Protein-like (EBPL)
Synonyms:	EBPL Antibody, ERP Antibody, Emopamil-binding protein-like Antibody, RP11-432M24.2 Antibody, EBRP Antibody, Emopamil binding protein-like Antibody
Host	EBPL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	EBPL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	EBPL antibody was raised against kLH-conjugated synthetic peptide (aa184-214) from the C-terminus of human EBPL.
Specificity:	Human EBPL
Epitope:	aa184-214
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences