

ESA / EPCAM Rabbit anti-Human Polyclonal (aa199-229) (FITC) Antibody - LS-C256746 - LSBio	
CatalogID:	LS-C256746
Target:	epithelial cell adhesion molecule (EPCAM)
Synonyms:	EPCAM Antibody, 323/A3 Antibody, ACSTD1 Antibody, 17-1A Antibody, CD326 Antibody, EGP Antibody, EGP34 Antibody, Epithelial glycoprotein Antibody, GA733 -2 Antibody, HNPCC8 Antibody, Ep-CAM Antibody, ESA Antibody, HEGP314 Antibody, KS 1/4 antigen Antibody, KS1/4 Antibody, KSA Antibody, M1S2 Antibody, MIC18 Antibody, MK-1 Antibody, MH99 Antibody, TROP1 Antibody, TACST-1 Antibody, TACSTD1 Antibody, CD326 antigen Antibody, CO-17A Antibody, DIAR5 Antibody, EGP-2 Antibody, EGP314 Antibody, EGP40 Antibody, Epithelial glycoprotein 314 Antibody, HEA125 Antibody, Ly74 Antibody, M4S1 Antibody, MOC31 Antibody
Host	EPCAM antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ESA / EPCAM antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ESA / EPCAM antibody was raised against kLH-conjugated synthetic peptide (aa199-229) from the C-terminus of human EpCAM.
Specificity:	Human ESA / EPCAM
Epitope:	aa199-229
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014 © 2014 LifeSpan BioSciences	