

ESR2 / ER Beta Rabbit anti-Human Polyclonal (aa451-480) (FITC) Antibody - LS-C256834 - LSBio	
CatalogID:	LS-C256834
Target:	estrogen receptor 2 (ER beta) (ESR2)
Synonyms:	ESR2 Antibody, ESTRB Antibody, ER Beta Antibody, ERbeta Antibody, ER-BETA Antibody, ESR-BETA Antibody, ESRB Antibody, Estrogen receptor 2 (ER beta) Antibody, NR3A2 Antibody, Estrogen receptor beta Antibody, Estrogen receptor beta 4 Antibody
Family / Subfamily:	NHR / NR3 Steroid receptor
Host	ESR2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ESR2 / ER Beta antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ESR2 / ER Beta antibody was raised against kLH-conjugated synthetic peptide (aa451-480) from the C-terminus of human ESR2.
Specificity:	Human ER Beta/Estrogen Receptor Beta
Epitope:	aa451-480
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014 © 2014 LifeSpan BioSciences	