

FGFR1 Rabbit anti-Human Polyclonal (pTyr766) (FITC) Antibody - LS-C257245 - LSBio

<b>CatalogID:</b>	LS-C257245
<b>Target:</b>	fibroblast growth factor receptor 1 (FGFR1)
<b>Synonyms:</b>	FGFR1 Antibody, BFGF receptor Antibody, BFGF-R Antibody, BFGF-R-1 Antibody, CD331 Antibody, Cek1 Antibody, CD331 antigen Antibody, FGFR1/PLAG1 fusion Antibody, FLT2 Antibody, Fgf-flg Antibody, FGFBR Antibody, FGFR-1 Antibody, Fms-like tyrosine kinase 2 Antibody, HBGFR Antibody, HH2 Antibody, Hydroxyaryl-protein kinase Antibody, FGF receptor 1 Antibody, Fms-related tyrosine kinase 2 Antibody, Proto-oncogene c-Fgr Antibody, BFGFR Antibody, CEK Antibody, FLT-2 Antibody, HBGF receptor Antibody, KAL2 Antibody, OGD Antibody
<b>Family / Subfamily:</b>	Protein Kinase / FGF Receptor
<b>Host</b>	FGFR1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	FGFR1 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	FGFR1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr766 of human FGFR.
<b>Specificity:</b>	Human FGFR1
<b>Epitope:</b>	pTyr766
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences