

FGFR1 Rabbit anti-Human Polyclonal (aa755-784) (FITC) Antibody - LS-C257257 - LSBio

CatalogID:	LS-C257257
Target:	fibroblast growth factor receptor 1 (FGFR1)
Synonyms:	FGFR1 Antibody, BFGF receptor Antibody, BFGF-R Antibody, BFGF-R-1 Antibody, CD331 Antibody, Cek1 Antibody, CD331 antigen Antibody, FGFR1/PLAG1 fusion Antibody, FLT2 Antibody, Fgf-flg Antibody, FGFBR Antibody, FGFR-1 Antibody, Fms-like tyrosine kinase 2 Antibody, HBGFR Antibody, HH2 Antibody, Hydroxyaryl-protein kinase Antibody, FGF receptor 1 Antibody, Fms-related tyrosine kinase 2 Antibody, Proto-oncogene c-Fgr Antibody, BFGFR Antibody, CEK Antibody, FLT-2 Antibody, HBGFR receptor Antibody, KAL2 Antibody, OGD Antibody
Family / Subfamily:	Protein Kinase / FGF Receptor
Host	FGFR1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	FGFR1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FGFR1 antibody was raised against kLH-conjugated synthetic peptide (aa755-784) from the C-terminus of Mouse Fgfr1.
Specificity:	Human FGFR1
Epitope:	aa755-784
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences