

GALNS / Chondroitinase Rabbit anti-Human Polyclonal (aa243-271) (PE) Antibody - LS-C257527 - LSBio

CatalogID:	LS-C257527
Target:	galactosamine (N-acetyl)-6-sulfate sulfatase (GALNS)
Synonyms:	GALNS Antibody, Chondroitinase Antibody, Chondroitinsulfatase Antibody, GALNAC6S Antibody, GalNAc6S sulfatase Antibody, GAS Antibody, MPS4A Antibody, Galactose-6-sulfate sulfatase Antibody
Host	GALNS antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	GALNS / Chondroitinase antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GALNS / Chondroitinase antibody was raised against kLH-conjugated synthetic peptide (aa243-271) from internal region of human GALNS.
Specificity:	Human GALNS / Chondroitinase
Epitope:	aa243-271
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences