

HCK Rabbit anti-Human Polyclonal (aa138-167) (FITC) Antibody - LS-C258329 - LSBio

CatalogID:	LS-C258329
Target:	hemopoietic cell kinase (HCK)
Synonyms:	HCK Antibody, Hematopoietic cell kinase Antibody, Hemopoietic cell kinase Antibody, JTK9 Antibody, Hck kinase Antibody, p59Hck Antibody, Tyrosine-protein kinase HCK Antibody, p59-hck Antibody, p61Hck Antibody, p60-hck Antibody
Family / Subfamily:	Protein Kinase / Src
Host	HCK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HCK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HCK antibody was raised against KLH-conjugated synthetic peptide (aa138-167) from the N-terminal region of human HCK.
Specificity:	Human HCK
Epitope:	aa138-167
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences