

HEYL Rabbit anti-Human Polyclonal (aa25-54) (FITC) Antibody - LS-C258439 - LSBio

CatalogID:	LS-C258439
Target:	hes-related family bHLH transcription factor with YRPW motif-like (HEYL)
Synonyms:	HEYL Antibody, HEY-like protein Antibody, HHRT3 Antibody, HRT-3 Antibody, HRT3 Antibody, HEY3 Antibody, BHLHb33 Antibody, HHeyL Antibody
Host	HEYL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HEYL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HEYL antibody was raised against kLH-conjugated synthetic peptide (aa25-54) from the N-terminal region of human HEYL.
Specificity:	Human HEYL
Epitope:	aa25-54
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences