

HOXB7 Rabbit anti-Rat Polyclonal (aa55-83) (PE) Antibody - LS-C258822 - LSBio

CatalogID:	LS-C258822
Target:	homeobox B7 (HOXB7)
Synonyms:	HOXB7 Antibody, Homeo box c1 protein Antibody, Homeobox protein HHO.C1 Antibody, Hox-2.3 Antibody, Homeo box 2C Antibody, Homeobox B7 Antibody, Homeobox protein Hox-B7 Antibody, HHO.C1 Antibody, Homeo box B7 Antibody, Homeobox protein Hox-2C Antibody, HOX2 Antibody, HOX2C Antibody
Host	HOXB7 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HOXB7 antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	HOXB7 antibody was raised against kLH-conjugated synthetic peptide (aa55-83) from internal region of human Rat Hoxb7.
Specificity:	Rat HOXB7.
Epitope:	aa55-83
Reactivity:	Rat, Rabbit
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences