

HOXD9 Rabbit anti-Rat Polyclonal (aa126-154) (PE) Antibody - LS-C258853 - LSBio

CatalogID:	LS-C258853
Target:	homeobox D9 (HOXD9)
Synonyms:	HOXD9 Antibody, Homeobox D9 Antibody, Homeobox protein Hox-5.2 Antibody, Hox-4.3 Antibody, Homeobox protein Hox-4C Antibody, Hox-4.3, mouse, homolog of Antibody, HOX4 Antibody, Homeo box D9 Antibody, Hox-5.2 Antibody, HOX4C Antibody, Homeobox protein Hox-D9 Antibody
Host	HOXD9 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HOXD9 antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	HOXD9 antibody was raised against kLH-conjugated synthetic peptide (aa126-154) from internal region of Rat Hoxd9.
Specificity:	Rat HOXD9
Epitope:	aa126-154
Reactivity:	Rat, Rabbit
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences