

DCC Rabbit anti-Mouse Polyclonal (aa642-670) (PE) Antibody - LS-C259167 - LSBio

CatalogID:	LS-C259167
Target:	deleted in colorectal carcinoma (DCC)
Synonyms:	DCC Antibody, CRC18 Antibody, Colorectal tumor suppressor Antibody, Colorectal cancer suppressor Antibody, IGDC1 Antibody, Netrin receptor DCC Antibody, Tumor suppressor protein DCC Antibody, CRCR1 Antibody
Host	DCC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	DCC antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	DCC antibody was raised against kLH-conjugated synthetic peptide (aa642-670) from internal region of human DCC.
Specificity:	Mouse DCC
Epitope:	aa642-670
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences