

Hamartin / TSC1 Rabbit anti-Human Polyclonal (Internal) (PE) Antibody - LS-C259664 - LSBio	
CatalogID:	LS-C259664
Target:	tuberous sclerosis 1 (TSC1)
Synonyms:	TSC1 Antibody, LAM Antibody, Tuberous sclerosis 1 protein Antibody, Tuberous sclerosis 1 Antibody, Hamartin Antibody, KIAA0243 Antibody, TSC Antibody, Tumor suppressor Antibody
Host	TSC1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Hamartin / TSC1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	Hamartin / TSC1 antibody was raised against kLH-conjugated synthetic peptide from internal region of human TSC1.
Specificity:	Human Hamartin / TSC1
Epitope:	Internal
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	