

IKBKB / IKK2 / IKK Beta Rabbit anti-Human Polyclonal (pSer682) (FITC) Antibody - LS-C259943 - LSBio	
<b>CatalogID:</b>	LS-C259943
<b>Target:</b>	inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase beta (IKBKB)
<b>Synonyms:</b>	IKBKB Antibody, Ikk 2 Antibody, IKK-B Antibody, IKK2 Antibody, IKKbeta Antibody, NFKB1KB Antibody, I-kappa-B kinase 2 Antibody, IKKB Antibody, I-kappa-B-kinase beta Antibody, Ikb kinase-beta Antibody, IKK Beta Antibody, IKK-beta Antibody
<b>Family / Subfamily:</b>	Protein Kinase / IKAPPAB
<b>Host</b>	IKBKB antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	IKBKB / IKK2 / IKK Beta antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	IKBKB / IKK2 / IKK Beta antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser682 of human IKKB.
<b>Specificity:</b>	Human IKBKB / IKK2 / IKK Beta
<b>Epitope:</b>	pSer682
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	