

IKBKB / IKK2 / IKK Beta Rabbit anti-Mouse Polyclonal (pSer701) (FITC) Antibody - LS-C260040 - LSBio	
CatalogID:	LS-C260040
Target:	inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase beta (IKBKB)
Synonyms:	IKBKB Antibody, Ikk 2 Antibody, IKK-B Antibody, IKK2 Antibody, IKKbeta Antibody, NFKB1KB Antibody, I-kappa-B kinase 2 Antibody, IKKB Antibody, I-kappa-B-kinase beta Antibody, Ikb kinase-beta Antibody, IKK Beta Antibody, IKK-beta Antibody
Family / Subfamily:	Protein Kinase / IKAPPAB
Host	IKBKB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IKBKB / IKK2 / IKK Beta antibody was raised against Mouse
Antigen Type:	Synthetic phospho-peptide
Immunogen:	IKBKB / IKK2 / IKK Beta antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser701 of mouse IKKB.
Specificity:	Mouse IKBKB / IKK2 / IKK Beta.
Epitope:	pSer701
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	