

**IL13 Rabbit anti-Human Polyclonal (aa117-146) (PE) Antibody - LS-C260084 - LSBio**

<b>CatalogID:</b>	LS-C260084
<b>Target:</b>	interleukin 13 (IL13)
<b>Synonyms:</b>	IL13 Antibody, Allergic rhinitis Antibody, ALRH Antibody, BHR1 Antibody, Interleukin-13 Antibody, NC30 Antibody, p600 Antibody, IL-13 Antibody, Interleukin 13 Antibody
<b>Family / Subfamily:</b>	Hormone / not assigned-Hormone
<b>Host</b>	IL13 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	IL13 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	IL13 antibody was raised against kLH-conjugated synthetic peptide (aa117-146) from the C-terminus of human IL13.
<b>Specificity:</b>	Human IL13 / IL-13
<b>Epitope:</b>	aa117-146
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences