

17B / IL-17B Rabbit anti-Human Polyclonal (aa46-74) (PE) Antibody - LS-C260103 - LSBio	
CatalogID:	LS-C260103
Target:	interleukin 17B (IL17B)
Synonyms:	IL17B Antibody, Cytokine Zcyto7 Antibody, Cytokine-like protein ZCYTO7 Antibody, IL-17B Antibody, Interleukin 17B Antibody, Interleukin 20 Antibody, Interleukin-17 beta Antibody, Interleukin-17B Antibody, NIRF Antibody, Interleukin-20 Antibody, ZCYTO7 Antibody
Family / Subfamily:	Interleukin
Host	IL17B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	17B / IL-17B antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	17B / IL-17B antibody was raised against kLH-conjugated synthetic peptide (aa46-74) from internal region of human IL17B.
Specificity:	Human 17B / IL-17B
Epitope:	aa46-74
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/29/2014 © 2014 LifeSpan BioSciences	