

IL4 Rabbit anti-Human Polyclonal (aa130-160) (PE) Antibody - LS-C260172 - LSBio

CatalogID:	LS-C260172
Target:	interleukin 4 (IL4)
Synonyms:	IL4 Antibody, B cell growth factor 1 Antibody, B-cell stimulatory factor 1 Antibody, B_cell stimulatory factor 1 Antibody, BCGF1 Antibody, BSF-1 Antibody, Binetrakin Antibody, BCGF-1 Antibody, Interleukin 4 Antibody, Interleukin-4 Antibody, Pitrakinra Antibody, BSF1 Antibody, IL-4 Antibody
Family / Subfamily:	Hormone / not assigned-Hormone
Host	IL4 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	IL4 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	IL4 antibody was raised against kLH-conjugated synthetic peptide (aa130-160) from the C-terminus of human IL4.
Specificity:	Human IL4 / IL-4
Epitope:	aa130-160
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences