

ILKAP Rabbit anti-Human Polyclonal (aa248-278) (PE) Antibody - LS-C260201 - LSBio

CatalogID:	LS-C260201
Target:	integrin-linked kinase-associated serine/threonine phosphatase (ILKAP)
Synonyms:	ILKAP Antibody, DKFZP434J2031 Antibody, ILKAP2 Antibody, ILKAP3 Antibody, pp2Cdelta Antibody, PP2C-DELTA Antibody
Family / Subfamily:	Protein Phosphatase / PP2C
Host	ILKAP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ILKAP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ILKAP antibody was raised against kLH-conjugated synthetic peptide (aa248-278) from internal region of human ILKAP.
Specificity:	Human ILKAP
Epitope:	aa248-278
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences